

# Mounting press

마운팅 프레스



**PRESS I**



**PRESS II**

The Press I specially developed for the production of mounted specimens.  
Top loading application of power and removal of mounted specimen.  
A safety locking system combined with precise screw locking clamp.  
All of operating displays on gauge and control entire cycle.  
One touch time and temperature setting button with autoamatic heating and cooling sequence system.  
Piston of cylinder moves up and down by hydraulic motor.

## PRESS I / II

MODEL	<b>Press I</b>	<b>Press II</b>
MOLD Q'TY	1	2
MOLD DIA.	(25.4, 31.7, 38.1, 40) mm	(25.4, 31.7, 38.1, 40) mm
OPERATION	LCD DISPLAY, BUTTON TOUCH	LCD DISPLAY, BUTTON TOUCH
ACCESSORY	POWDER 500g, SPOON	POWDER 500g, SPOON
TYPE	FULLY AUTOMATIC CYCLE, HYDRAULIC CYLINDER	FULLY AUTOMATIC CYCLE, HYDRAULIC CYLINDER
POWER	220V, 1PH, 1 KW	220V, 1PH, 1 KW